



## PERFECTO CLASSIC MODEL HEARING AID

ASUK-9800-1000

- Easy to use
- Good sound quality
- Adjustable volume control
- Designed for use when using the telephone
- Lightweight to wear
- Protective storage case for when not in use
- Easy to follow instructions in Arabic and English
- 12 month Warranty
- Full after-sales service support

### TECHINICAL SPECIFICATION

MAX SSPL90	128dB SPL
HF Average SSPL90	120dB SPL
Peak Gain	54dB
HFA Full On Gain	48dB
Frequency Range	200Hz – 5300Hz
THD at	500Hz: 3.8%
	800Hz: 1.2%
	1600Hz: 0.2%
Eq. Input Noise Level	28dB
Battery Current Drain	1.4mA
HFA Splits	93dB

