

# Streamline2000Systems



### STREAMLINE 2000 SYSTEM

The Streamline 2000 System is the way forward for diagnostic instruments. Each option has been designed to deal with the everyday patient ophthalmic and aural examinations carried out by trained medical staff. Early diagnosis of a problem can identify the appropriate treatment required – which in many cases leads to a complete recovery.

With over 27 variations and combinations of instruments available in this range there is an option to suit almost every medical treatment centre, from the large hospital to the emergency crews responding to distress calls.

Pages 4 to 12 show the many options that are available with the Streamline 2000 System and there is more product information shown in the Specification Section. – see page 13 & 15



## CONTENTS

Description	Page
Mini Fibre Optic Otoscope Set	4
3.5V Fibre Optic Desk Top Otoscope Set	4
3.5V Fibre Optic Rechargeable Otoscope Set	4
3.5V Fibre Optic Wall Mounted Otoscope Set	4
Direct Mini Ophthalmoscope Set with 4 Apertures	5
Direct Mini Ophthalmoscope Set with 5 Apertures	5
3.5V Coaxial Ophthalmoscope Desk Top Set	5
3.5V Coaxial Ophthalmic Rechargeable Set	5
3.5V Coaxial Ophthalmoscope Wall Mounted Set	6
3.5V Streak Retinoscope Desk Top Set	6
3.5V Streak Retinoscope Rechargeable Set	6
3.5V Streak Retinoscope Wall Mounted Set	6
GP Combi Otoscope & Ophthalmoscope Desk Top Set	7
GP Combi Otoscope & Ophthalmoscope Rechargeable Set	7
GP Combi Otoscope & Ophthalmoscope Wall Mounted Set	7
Ophthalmoscope & Retinoscope Combi Desk Top	7
Ophthalmoscope & Retinoscope Combi Rechargeable Set	8
Ophthalmoscope & Retinoscope Combi Wall Mounted Set	8
GP Twin Handle Otoscope & Ophthalmoscope Rechargeable Set	8
Ophthalmoscope & Retinoscope Twin Handle Rechargeable Set	8
Ophthalmoscope & Retinoscope Twin Handle Wall Mounted Set	9
GP Basic Twin Handle Otoscope & Ophthalmoscope Wall Mounted Unit	9
GP Basic Plus Twin Handle Otoscope with Speculum Dispenser & Ophthalmoscope Wall Mounted Unit	9
GP Basic Extra Plus Twin Handle Otoscope with Speculum Dispenser, Ophthalmoscope & Sphygmomanometer Wall Mounted Unit	9
GP Pro Plus Twin Handle Otoscope with Speculum Dispenser, Ophthalmoscope, Ear Thermometer & Sphygmomanometer Wall Mounted Unit	10
Mini Diagnostic Set with 4 Apertures	10
Mini Diagnostic Set with 5 Apertures	10
GP Basic 3.5V Fibre Optic Otoscope Set with Hard Case	10
Hard Case Ophthalmoscope Set	11
GP Hard Case Basic Otoscope Set// Plus Throat Illuminator	11
Otoscope Hard Case Illuminator Set	11
GP Hard Case Combi Rechargeable Diagnostic Set	11
GP Hard Case Combi Rechargeable Diagnostic Plus Throat Illuminator Set	12
GP Hard Case Combi Rechargeable Procedure Set	12
Hard Case Ophthalmic - Retinoscope Set	12
ENT Diagnostic Basic Instrument Set	12

## Code Description & Specification

2000-1001 Mini Fibre Optic Set	13
2000-1005 Direct Mini Ophthalmoscope 4 Aperture Set	13
2000-1006 Direct Mini Ophthalmoscope 5 Aperture Set	13
2000-2001 3.5V Fibre Optic Otoscope Head	13
2000-2003 3.5V Coaxial Ophthalmoscope Head	13
2000-2004 3.5V Retinoscope Head	13
2000-4002 Single Charger Unit	14
2000-4004 Twin Charger Unit	14
2000-4007 Rechargeable Handle	14
2000-4008 Desk Top Power Unit	14
2000-4009 Single Wall Mount Transformer	14
2000-4010 Twin Wall Mount Transformer	14
2000-6501 Throat Illuminator	15
2000-6502 Fibre Optic Nose & Throat Illuminator	15
2000-6503 Fibre Optic Illuminated Nasal Speculum	15
2000-6504 Fibre Optic Illuminated Bi-Valve Nasal Speculum	15
2000-6505 Fibre Optic Illuminated Laryngeal Mirror	15
2000-6506 Fibre Optic Illuminated Operating Otoscope	15

2000-1001

### Streamline Mini Fibre Optic Otoscope Set

- 2.5 x Magnification Swivel wide lens to allow for minor procedures.
- Insufflation Port.
- Supplied complete with starter pack of 10 re-usable speculum ( 5 pieces of 2.5mm & 4.0mm).
- Hard Case.
- Handle accepts AA or Ni-MH batteries



2000-1002

### Streamline Desktop Fibre Optic Otoscope Set

- 3.5v Fiber Optic Otoscope - Halogen Lamp.
- 2 x Magnification Swivel wide angle lens to allow for minor procedures.
- Insufflation Port.
- Supplied complete with starter pack of 10 re-usable speculum ( 5 pieces of 2.7mm & 4.2mm).
- Desktop Power Unit with Handle and Power Adapter plus wall plug.



2000-1003

### Streamline Rechargeable Fibre Optic Otoscope Set

- 3.5v Fiber Optic Otoscope - Halogen Lamp.
- 2 x Magnification Swivel wide angle lens to allow for minor procedures.
- Insufflation Port.
- Supplied complete with starter pack of 10 re-usable speculum ( 5 pieces of 2.7mm & 4.2mm).
- Rechargeable Handle with battery.
- Charger Unit with Power Adaptor plus wall plug



2000-1004

### Streamline Wall Mounted Fibre Optic Otoscope Set

- 3.5v Fiber Optic Otoscope
- Halogen Lamp.
- 2 x Magnification Swivel wide angle lens to allow for minor procedures.
- Insufflation Port.
- Supplied complete with starter pack of 10 re-usable speculum ( 5 pieces of 2.7mm & 4.2mm).
- Handle, Single Wall Mount Transformer with power cable plus plug





2000-1005

### Streamline Direct Mini Ophthalmoscope Set

4 Standard Apertures

Red-free filter, hemislot, small spot, large spot,

- 19 Lenses, 4 Apertures.
- Hard Case.
- Handle accepts AA or Ni-MH batteries.



2000-1006

### Streamline Direct Mini Ophthalmoscope Set

5 Standard Apertures

Small spot, large spot, red-free filter, Slit, Fixation aperture.

- 19 Lenses, 5 Apertures.
- Hard Case.
- Handle accepts AA or Ni-MH batteries.



2000-1007

### Streamline Desktop Ophthalmoscope Set

- 3.5v Coaxial Ophthalmoscope .
- With 28 lenses and 18 aperture/filter combinations.
- Desktop Power Unit with Handle and Power Adapter plus wall plug.



2000-1008

### Streamline Rechargeable Ophthalmoscope Set

- 3.5v Coaxial Ophthalmoscope.
- With 28 lenses and 18 aperture/filter combinations.
- Rechargeable Handle with battery.
- Charger Unit with Power Adaptor plus wall plug.



2000-1009

### **Streamline Wall Mounted Ophthalmoscope Set**

- 3.5v Coaxial Ophthalmoscope.
- With 28 lenses and 18 aperture/filter combinations.
- Handle.
- Single Wall Mount Transformer with power cable plus plug.



2000-1010

### **Streamline Desktop Streak Retinoscope Set**

- 3.5v Streak Retinoscope.
- Desktop Power Unit with Handle.
- Power Adapter plus wall plug.



2000-1011

### **Streamline Rechargeable Streak Retinoscope Set**

- 3.5v Streak Retinoscope.
- Rechargeable Handle with battery.
- Charger Unit with Power Adaptor plus wall plug



2000-1012

### **Streamline Wall Mounted Streak Retinoscope Set**

- 3.5v Streak Retinoscope.
- Handle.
- Single Wall Mount Transformer with power cable plus plug.



2000-10015

### Streamline GP Combi Desktop Diagnostic Set

- 3.5v Fiber Optic Otoscope
- Halogen Lamp.
- 2 x Magnification Swivel wide angle lens to allow for minor procedures.
- Insufflation Port.
- Supplied complete with starter pack of 10 re-usable speculum ( 5 pieces of 2.7mm & 4.2mm).
- 3.5v Coaxial Ophthalmoscope with 28 lenses and 18 aperture/filter combinations.
- Desktop Power Unit with Handle and Power Adapter plus wall plug



2000-1016

### Streamline GP Combi Rechargeable Diagnostic Set

- 3.5v Fiber Optic Otoscope.
- Halogen Lamp.
- 2 x Magnification Swivel wide angle lens to allow for minor procedures.
- Insufflation Port.
- Supplied complete with starter pack of 10 re-usable speculum ( 5 pieces of 2.7mm & 4.2mm).
- 3.5v Coaxial Ophthalmoscope with 28 lenses and 18 aperture/filter combinations.
- Rechargeable Handle with battery.
- Single Charger Unit with Power Adaptor plus wall plug.



2000-1017

### Streamline GP Combi Wall Mounted Diagnostic Set

- 3.5v Fiber Optic Otoscope.
- Halogen Lamp.
- 2 x Magnification Swivel wide angle lens to allow for minor procedures.
- Insufflation Port.
- Supplied complete with starter pack of 10 re-usable speculum ( 5 pieces of 2.7mm & 4.2mm).
- 3.5v Coaxial Ophthalmoscope with 28 lenses and 18 aperture/filter combinations.
- Handle.
- Single Wall Mount Transformer with power cable plus plug.



2000-1018

### Streamline Ophthalmic Combi Desktop Diagnostic Set

- 3.5v Coaxial Ophthalmoscope
- With 28 lenses and 18 aperture/filter combinations.
- 3.5v Streak Retinoscope.
- Desktop Power Unit with Handle and Power Adapter plus wall plug



2000-1019

### Streamline Ophthalmic Combi Rechargeable Diagnostic Set

- 3.5v Coaxial Ophthalmoscope.
- With 28 lenses and 18 aperture/filter combinations.
- 3.5v Streak Retinoscope. Rechargeable Handle with battery.
- Single Charger Unit with Power Adaptor plus wall plug.



2000-1020

### Streamline Ophthalmic Combi Wall Mounted Diagnostic Set

- 3.5v Coaxial Ophthalmoscope.
- With 28 lenses and 18 aperture/filter combinations.
- 3.5v Streak Retinoscope. Handle.
- Single Wall Mount Transformer with power cable plus plug.



2000-1021

### Streamline GP Twin Base Rechargeable Diagnostic Set

- 3.5v Fiber Optic Otoscope
- Halogen Lamp.
- 2 x Magnification Swivel wide angle lens to allow for minor procedures.
- Insufflation Port.
- Supplied complete with starter pack of 10 re-usable speculum ( 5 pieces of 2.7mm & 4.2mm).
- 3.5v Coaxial Ophthalmoscope
- With 28 lenses and 18 aperture/filter combinations.
- 2 Rechargeable Handles with batteries.
- Twin Charger Unit with Power Adaptor plus wall plug



2000-1022

### Streamline Ophthalmic Twin Base Rechargeable Diagnostic Set

- 3.5v Coaxial Ophthalmoscope
- with 28 lenses and 18 aperture/filter combinations.
- 3.5v Streak Retinoscope.
- 2 Rechargeable Handles with batteries.
- Twin Charger Unit with Power Adaptor plus wall plug





2000-1023

### Streamline Ophthalmic Twin Wall Mounted Diagnostic Set

- 3.5v Coaxial Ophthalmoscope
- with 28 lenses and 18 aperture/filter combinations.
- 3.5v Streak Retinoscope.
- 2 Handles,
- Twin Wall Mount Transformer with power cable plus plug.



2000-1024

### Streamline GP Basic Twin Wall Mounted Diagnostic Set

- 3.5v Fiber Optic Otoscope
- Halogen Lamp.
- 2 x Magnification Swivel wide angle lens to allow for minor procedures.
- Insufflation Port.
- Supplied complete with starter pack of 10 re-usable speculum ( 5 pieces of 2.7mm & 4.2mm).
- 3.5v Coaxial Ophthalmoscope
- with 28 lenses and 18 aperture/filter combinations.
- 2 Handles, Twin Wall Mount Transformer
- with power cable plus plug.



2000-1025

### Streamline GP Basic Plus Twin Wall Mounted Diagnostic Set

- 3.5v Fiber Optic Otoscope
- Halogen Lamp.
- 2 x Magnification Swivel wide angle lens to allow for minor procedures.
- Insufflation Port.
- Supplied complete with starter pack of 10 re-usable speculum ( 5 pieces of 2.7mm & 4.2mm).
- Wall Mounted Speculum Dispenser.
- 3.5v Coaxial Ophthalmoscope with 28 lenses and 18 aperture/filter combinations.
- 2 Handles, Twin Wall Mount Transformer with power cable plus plug



2000-10026

### Streamline GP Basic Extra Plus Twin Wall Mounted Diagnostic Set

- 3.5v Fiber Optic Otoscope
- Halogen Lamp.
- 2 x Magnification Swivel wide angle lens to allow for minor procedures.
- Insufflation Port.
- Supplied complete with starter pack of 10 re-usable speculum ( 5 pieces of 2.7mm & 4.2mm).
- Wall Mounted Speculum Dispenser.
- 3.5v Coaxial Ophthalmoscope with 28 lenses and 18 aperture/filter combinations.
- Sphygmomanometer.
- 2 Handles, Twin Wall Mount Transformer with power cable plus plug



2000-1027

## Streamline GP Pro Plus Twin Wall Mounted Diagnostic Set

- 3.5v Fiber Optic Otoscope
- Halogen Lamp.
- 2 x Magnification Swivel wide angle lens to allow for minor procedures.
- Insufflation Port.
- Supplied complete with starter pack of 10 re-usable speculum ( 5 pieces of 2.7mm & 4.2mm).
- Wall Mounted Speculum Dispenser.
- 3.5v Coaxial Ophthalmoscope with 28 lenses and 18 aperture/filter combinations.
- Sphygmomanometer.
- Ear Thermometer.
- 2 Handles, Twin Wall Mount Transformer with power cable plus plug



2000-5001

## Streamline Mini Diagnostic Hard Case Set

4 Standard Apertures



Red-free filter, hemispherical spot, small spot, large spot,

- Fibre Optic Otoscope.
- 2.5 x Magnification Swivel wide angle lens to allow for minor procedures.
- Insufflation Port.
- Supplied complete starter pack of 10 re-usable speculum ( 5 pieces of 2.5mm & 4.0mm).
- Ophthalmoscope 19 Lenses,
- 4 Apertures.
- Hard Case.
- Handle accepts AA or Ni-MH batteries



2000-5002

## Streamline Mini Diagnostic Hard Case Set

5 Standard Apertures



Small spot, large spot, red-free filter, Slit, Fixation aperture.

- Fibre Optic Otoscope.
- 2.5 x Magnification Swivel wide angle lens to allow for minor procedures.
- Insufflation Port.
- Supplied complete starter pack of 10 re-usable speculum ( 5 pieces of 2.5mm & 4.0mm).
- Ophthalmoscope 19 Lenses,
- 5 Apertures.
- Hard Case.
- Handle accepts AA or Ni-MH batteries



2000-5003

## Streamline GP Hard Case Basic Otoscope Set

- 3.5v Fiber Optic Otoscope –
- Halogen Lamp.
- 2 x Magnification Swivel wide angle lens to allow for minor procedures.
- Insufflation Port.
- Supplied complete with starter pack of 10 re-usable speculum ( 5 pieces of 2.7mm & 4.2mm).
- Rechargeable Handle with battery.
- Hard Case .
- Charger Unit with Power Adaptor plus wall plug



2000-5004

### Streamline Hard Case Ophthalmoscope Set

- 3.5v Coaxial Ophthalmoscope with 28 lenses and 18 aperture/filter combinations.
- Rechargeable Handle with battery.
- Hard Case.
- Charger Unit with Power Adaptor plus wall plug



2000-5005

### Streamline GP Hard Case Basic Otoscope Set Plus Throat Illuminator

- 3.5v Fiber Optic Otoscope
- Halogen Lamp.
- 2 x Magnification Swivel wide angle lens to allow for minor procedures.
- Insufflation Port.
- Throat Illuminator
- Supplied complete with starter pack of 10 re-usable speculum ( 5 pieces of 2.7mm & 4.2mm).
- Rechargeable Handle with battery.
- Hard Case .
- Charger Unit with Power Adaptor plus wall plug



2000-5006

### Streamline Otoscope Hard Case Illuminator Set

- 3.5v Fiber Optic Otoscope
- Halogen Lamp.
- 2 x Magnification Swivel wide angle lens to allow for minor procedures.
- Insufflation Port.
- Supplied complete with starter pack of 10 re-usable speculum ( 5 pieces of 2.7mm & 4.2mm).
- Throat Illuminator. Fibre Optic Illuminator. Fibre Optic Illuminated Nasal Speculum.
- Fibre Optic Illuminated Bi-Valve Nasal Speculum.
- Fibre Optic Laryngeal Mirror.
- Rechargeable Handle with battery.
- Hard Case .Charger Unit with Power Adaptor plus wall plug



2000-5007

### Streamline GP Hard Case Combi Rechargeable Diagnostic Set

- 3.5v Fiber Optic Otoscope
- Halogen Lamp.
- 2 x Magnification Swivel wide angle lens to allow for minor procedures.
- Insufflation Port.
- Supplied complete with starter pack of 10 re-usable speculum ( 5 pieces of 2.7mm & 4.2mm).
- 3.5v Coaxial Ophthalmoscope with 28 lenses and 18 aperture/filter combinations.
- Rechargeable Handle with battery.
- Single Charger Unit with Power Adaptor plus wall plug



2000-5008

### **Streamline GP Hard Case Combi Rechargeable Diagnostic Plus Throat Illuminator Set**

- 3.5v Fiber Optic Otoscope
- Halogen Lamp.
- 2 x Magnification Swivel wide angle lens to allow for minor procedures.
- Insufflation Port.
- Supplied complete with starter pack of 10 re-usable speculum ( 5 pieces of 2.7mm & 4.2mm).
- 3.5v Coaxial Ophthalmoscope with 28 lenses and 18 aperture/filter combinations.
- Throat Illuminator.
- Rechargeable Handle with battery. Single Charger Unit with Power Adaptor plus wall plug



2000-5009

### **Streamline GP Hard Case Combi Rechargeable Procedure Set**

- 3.5v Fibre Optic Procedure Otoscope with 5 Reusable Speculum 2mm/3mm/4mm/5mm.9mm
- Halogen Lamp.
- 2 x Magnification Swivel lens to allow for minor procedures.
- Insufflation Port.
- 3.5v Coaxial Ophthalmoscope with 28 lenses and 18 aperture/filter combinations.
- Throat Illuminator.
- Rechargeable Handle with battery.
- Single Charger Unit with Power Adaptor plus wall plug.



2000-5010

### **Streamline Hard Case Ophthalmic – Retinoscope Set**

- 3.5v Coaxial Ophthalmoscope with 28 lenses and 18 aperture/filter combinations.
- 3.5v Streak Retinoscope.
- Throat Illuminator.
- Rechargeable Handle with battery.
- Single Charger Unit with Power Adaptor plus wall plug



6001-1040

### **ENT Diagnostic Basic Instrument Set**

- Black Light Stem LED Cold Lamp
- Tongue Depressor – Stainless Steel
- Holder for Wooden Tongue Depressor
- Laryngeal Mirror No3 (20mm diameter) with Light Stem Fitting
- Laryngeal Mirror No4 (22mm diameter) with Light Stem Fitting
- Vinyl Zipper Case
- Fits all 3.5v Handles



## Code Description & Specification



**2000-1001 Streamline Mini Fibre Optic Set** – see page 4  
**2000-7001 Reusable Speculum 2.5mm**  
**2000-7002 Reusable Speculum 4.0mm**  
**2000-3001 2.5V Replacement Lamp**



**2000-1005 Streamline Direct Mini Ophthalmoscope Set** – see page 5  
**2000-1006 Streamline Direct Mini Ophthalmoscope Set** – see page 5  
**2000-3004 2.5V Replacement Lamp**



**2000-2001 3.5V Fibre Optic Otoscope Head**  
 2 x Magnification  
 Sealed System for Pneumatic Otoscopy  
 Wide Angle Swivel Lens to suit procedures  
 Interchangeable with all 3.5V Handles  
**2000-7003 Reusable Speculum 2.7mm**  
**2000-7004 Reusable Speculum 4.2mm**  
**7001-4005 Disposable Speculum 2.75mm**  
**7001-4006 Disposable Speculum 4.2mm**  
**2000-3002 3.5V Replacement Lamp**



**2000-2003 3.5V Coaxial Direct Ophthalmoscope Head**  
 Coaxial Optics allow easy entry into undilated pupils  
 Enables a more accurate diagnosis  
 3.5V Halogen Lamp provides quality light output for true tissue colour and long lasting consistent illumination  
 28 Lenses from: 0D, +/-1D, +/-2D, +/-3D, +/-4D, +/-5D, +/-6D, +/-7D, +/-8D, +/-9D, +/-10D, +12D, +/-15D, +/-20D, -25D, +40D  
 18 Unique Aperture/Filter Combinations giving greater versatility –add colourways as shown on page 5 of pdf brochure  
 Interchangeable with all 3.5V Handles  
**2000-3004 3.5V Replacement Lamp**



**2000-2004 3.5V Streak Retinoscope Head**  
 Fast, accurate and easy to use  
 Excellent for determining the axis of astigmatism  
 Superior streak qualities for easier and quicker refractions, even in patients with small pupils  
 Bright lamp ensures easy recognition of neutralism  
 Multi-coated optics  
 Single control for vergence and rotation  
 Lowest position gives a divergent beam  
 Dust proof and maintenance free  
 Metal bayonet connector designed to give Retinoscope a longer working operation  
 Interchangeable with all 3.5V Handles  
**2000-3002 3.5V Replacement Lamp**



## Code Description & Specification



### **2000-4002 Single Charger Unit**

### **2000-4005 Optional Single Charger Unit Wall Mounting Bracket**

### **2000-4004 Twin Charger Unit**

### **2000-4006 Optional Twin Charger Unit Wall Mounting Bracket**

Compact portable units providing convenient recharging  
 Switch Mode: Single-Cell/Two Cell Li-ion/Li Polymer Battery Charger.  
 Precharge- CC-CV Operation model  
 Short circuit, reverse voltage and over voltage protection  
 Low Battery Status Indicator by Yellow-Green LED  
 Supplied complete with Power Adaptor to suit local requirements



### **2000-4007 Rechargeable Handle**

Half the weight of traditional handles  
 Comfortable fit in the hand avoids strain  
 Continuous supply of adjustable brightness  
 Made from high-impact ABS material  
 Complete with rechargeable battery  
 Fits all 3.5V instrument heads



### **2000-4008 Desk Top Power Unit**

Provides a constant 3.5V output voltage  
 Handle enables continuous adjustable brightness  
 Designed to give a comfortable and lightweight fit in the hand  
 Fits all 3.5V instrument heads  
 Automatic On/Off mechanism when the handle is removed or replaced from the unit  
 Supplied complete with Power Cable to suit all local requirements



### **2000-4009 Single Wall Mount Transformer**

### **2000-4010 Twin Wall Mount Transformer**

Provides a constant 3.5V output independent of voltage fluctuations  
 Automatic On/Off switching when the handle is either removed or replaced from its holder saving energy and increasing the working life of the lamps  
 Continuous brightness control on the handle  
 Designed to give a comfortable and lightweight fit in the hand  
 Quality spiral cords giving a comfortable working radius of upto 3 metres  
 Handles fit all 3.5V instrument heads  
 Supplied complete with Power Cable to suit all local requirements

## Code Description & Specification



2000-6501

### Throat Illuminator

Halogen Lamp. Fits all 3.5V Handles



2000-6502

### Fibre Optic Nose & Throat Illuminator

Cold white light with no risk of heat burns to patient.  
Fits all 3.5V Handles.



2000-6503

### Fibre Optic Illuminated Nasal Speculum

Cold white light with no risk of heat burns to patient.  
Fits all 3.5V Handles



2000-6504

### Fibre Optic Illuminated Bi-Valve Nasal Speculum.

Cold white light with no risk of heat burns to patient.  
Fits all 3.5V Handles



2000-6505

### Fibre Optic Illuminated Laryngeal Mirror

Cold white light with no risk of heat burns to patient.  
Fits all 3.5V Handles



2000-6506

### Fibre Optic Illuminated Operating Otoscope

Supplied with 5 Sizes of reusable speculum – 2mm/3mm/4mm/5mm/9mm  
Swivel magnifying lens to give added increased space for when instruments used on  
minor procedures.  
Fits all 3.5V Handles



Telephone: +44 (0) 1202 579 449  
Fax: +44 (0) 1202 552 444  
Email: [info@avondalesurgicaluk.com](mailto:info@avondalesurgicaluk.com)  
[www.avondalesurgicaluk.com](http://www.avondalesurgicaluk.com)

Unit A, Building 402 Target Road,  
Aviation Business Park,  
Bournemouth Airport Christchurch.  
Dorset BH23 6NW UK



All Medical Devices are CE marked in  
accordance with the Medical Device  
Directive (93/42/EEC).